

YOUR INDIVIDUAL OPTIONS

INTEGO CS

Options	INTEGO CS
Patient chair	
2 treatment positions, S and Q position	●
Number of user-specific settings	2
Central junction box	●
Chair lift system 140 kg	■
Chair lift system 185 kg	□
Flat headrest	■
Double articulating headrest	□
Premium upholstery	■
4-way foot switch	○
Booster cushion for small patients	○
Armrest, right	○
Armrest, left	○
Valve set for external devices	○
Pressure reduction kit	○
Cover installation kit	○
Separate junction box	○
Base plate for uneven floors, adapts for C-series installation to be replaced	○
Dentist element	
Height-adjustable support arm with pneumatic brake	●
5 instrument positions	●
Hydro-pneumatic supply for 2 drill drives w/o fiber optics	●
Pneumatic instrument foot control	■
C ⁺ electronic foot control	□ A, I
EasyPad user interface	●
EasyTouch user interface	–
SPRAYVIT E with light	□ A, I
SPRAYVIT E with light (water on the left)	□ A, I
Standard 3-way syringe	□ S
Turbine hose with light	○ 2xS, A, I
Borden turbine hose	○
BL ISO E micro-motor with hose	○ A, I
BL ISO E micro-motor with hose (ApexLocator ready)	–
Second BL ISO E micro-motor (ApexLocator ready)	–
SIROSONIC L scaler	○ I
Third drill drive	○
Satelec Mini LED curing light	–
Intra-oral light source kit for turbine and SPRAYVIT E	○ S, A, I
Torque control kit	○
ApexLocator	–
ENDO 6:1 contra-angle handpiece	○
Intraoral X-ray viewer	○
Purge and manual sanitization function	○
Purge function	–
Additional CS tray holder for 1 standard size tray	○
Additional CS tray holder for 2 standard size trays	○
Silicone pad for CS dentist element	○
Removable silicone sleeves for handles	○

Options	INTEGO CS
Assistant element	
Compact version, swivel mounted directly on the water unit	■
Comfort version, on a support arm with wide swivel range	□
4 handpiece positions	●
Saliva ejector and HVE	●
SPRAYVIT E with light	□
SPRAYVIT E with light (water on the left)	□
Standard 3-way syringe	□
2nd HVE kit	○
Satelec Mini LED curing light	○
Touch pad	○
Light source for SPRAYVIT E	○
Silicone pad for Comfort assistant element	○
Water unit	
Compact design, prepared for fresh water bottle	●
Comfort design, prepared for integrated disinfection device	–
Cuspidor with manual swivel function	●
Reduced price without cuspidor	□
Prepared for standard wet suction	■
Prepared for Standard wet suction w/automatic suction control	□
Multi-station valve	□
Dürr cuspidor valve	□
Automatic separator	□
Amalgam separator	□
Air jet pump (Air Venturi)	□
Fresh water bottle w/o switch to mains water supply	□
Fresh water bottle with switch to mains water supply	□
Integrated disinfection system with permanent waterline disinfection and automatic sanitization function	–
Water heater	○
Water flush suction hose cleaning system	○
Chemical suction hose cleaning system	–
Operating light options	
LEDlight	□ S, A
LEDview	□ I
Adaptation for HELIODENT Plus unit model	○
Patient communication	
USB interface in the dentist or assistant element	○ SiR, SiS, SiA
SIROCAM F	○ SiS
SIROCAM AF	○ SiA
Monitor mount on the light support column	○ SiS, SiA
22" AC monitor (with integrated power supply unit)	○ SiS*, SiA
Network interface	○ SiR, SiS, SiA

- Included in the basic unit
- Option
- Variation, part of standard unit
- Optional variation
- Not possible/applicable

- SiR SIVISION ready package
- SiS SIVISION Start package with/without monitor
- SiS* SIVISION Start package with monitor
- SiA SIVISION Advance package
- S START/START pro package
- A ADVANCE/ADVANCE pro package
- I INNOVATION/INNOVATION pro package
- E ENDO pro package