

YOUR INDIVIDUAL OPTIONS

INTEGO pro CS

| Options | INTEGO Pro CS |
|---|---------------|
| Patient chair | |
| 2 treatment positions, S and 0 position | ● |
| Number of user-specific settings | 4 |
| Central junction box | ● |
| Chair lift system 140 kg | – |
| Chair lift system 185 kg | ● |
| Flat headrest | ■ |
| Double articulating headrest | □ |
| Premium upholstery | ■ |
| 4-way foot switch | ● |
| Booster cushion for small patients | ○ |
| Armrest, right | ○ |
| Armrest, left | ○ |
| Valve set for external devices | ○ |
| Pressure reduction kit | ○ |
| Cover installation kit | ○ |
| Separate junction box | ○ |
| Base plate for uneven floors, adapts for C-series installation to be replaced | ○ |
| Dentist element | |
| Height-adjustable support arm with pneumatic brake | ● |
| 5 instrument positions | ● |
| Hydro-pneumatic supply for 2 drill drives w/o fiber optics | ● |
| Pneumatic instrument foot control | ■ |
| C* electronic foot control | □ S, A, I |
| EasyPad user interface | ■ |
| EasyTouch user interface | □ A, I |
| SPRAYVIT E with light | □ S, A, I |
| SPRAYVIT E with light (water on the left) | □ S, A, I |
| Standard 3-way syringe | – |
| Turbine hose with light | ○ S, A, I |
| Borden turbine hose | – |
| BL ISO E micro-motor with hose | – |
| BL ISO E micro-motor with hose (ApexLocator ready) | ○ S, A, I |
| Second BL ISO E micro-motor (ApexLocator ready) | ○ I |
| SIROSONIC L scaler | ○ A, I |
| Third drill drive | ○ I |
| Satelec Mini LED curing light | ○ |
| Intra-oral light source kit for turbine and SPRAYVIT E | ● |
| Torque control kit | ○ E |
| ApexLocator | ○ E |
| ENDO 6:1 contra-angle handpiece | ○ E |
| Intraoral X-ray viewer | ○ |
| Purge and manual sanitization function | – |
| Purge function | ● |
| Additional CS tray holder for 1 standard size tray | ○ |
| Additional CS tray holder for 2 standard size trays | ○ |
| Silicone pad for CS dentist element | ○ |

| Options | INTEGO Pro CS |
|--|-----------------|
| Removable silicone sleeves for handles | ○ |
| Assistant element | |
| Compact version, swivel mounted directly on the water unit | – |
| Comfort version, on a support arm with wide swivel range | ● |
| 4 handpiece positions | ● |
| Saliva ejector and HVE | ● |
| SPRAYVIT E with light | □ I |
| SPRAYVIT E with light (water on the left) | □ I |
| Standard 3-way syringe | – |
| 2nd HVE kit | ○ |
| Satelec Mini LED curing light | ○ |
| Touch pad | ○ |
| Light source for SPRAYVIT E | ● |
| Silicone pad for Comfort assistant element | ○ |
| Water unit | |
| Compact design, prepared for fresh water bottle | – |
| Comfort design, prepared for integrated disinfection device | ● |
| Cuspidor with manual swivel function | ● |
| Reduced price without cuspidor | – |
| Prepared for standard wet suction | – |
| Prepared for Standard wet suction w/automatic suction control | – |
| Multi-station valve | ■ |
| Dürr cuspidor valve | □ |
| Automatic separator | □ |
| Amalgam separator | □ |
| Air jet pump (Air Venturi) | – |
| Fresh water bottle w/o switch to mains water supply | – |
| Fresh water bottle with switch to mains water supply | – |
| Integrated disinfection system with permanent waterline disinfection and automatic sanitization function | ○ |
| Water heater | ● |
| Water flush suction hose cleaning system | ■ |
| Chemical suction hose cleaning system | □ |
| Operating light options | |
| LEDlight | □ S |
| LEDview | □ A, I |
| Adaptation for HELIODENT Plus unit model | ○ |
| Patient communication | |
| USB interface in the dentist or assistant element | ○ SiR, SiS, SiA |
| SIROCAM F | ○ SiS |
| SIROCAM AF | ○ SiA |
| Monitor mount on the light support column | ○ SiS, SiA |
| 22" AC monitor (with integrated power supply unit) | ○ SiS*, SiA |
| Network interface | ○ SiR, SiS, SiA |

- Included in the basic unit
- Option
- Variation, part of standard unit
- Optional variation
- Not possible/applicable

- SiR SIVISION ready package
- SiS SIVISION Start package with/without monitor
- SiS* SIVISION Start package with monitor
- SiA SIVISION Advance package
- S START/START pro package
- A ADVANCE/ADVANCE pro package
- I INNOVATION/INNOVATION pro package
- E ENDO pro package